

# Action 3 Jnr

CREATIVE THINKING



Invacare's Action<sup>3</sup> Junior is a lightweight, foldable paediatric wheelchair designed for children aged between 3 and 15 years. Action<sup>3</sup> Junior has been developed to match the individual needs of the child and can grow as they grow. The Action<sup>3</sup> Junior comes in attractive, modern styles, designed to appeal to kids.



Yes, you can.®

Making Life's  
Experiences Possible



# Action 3 Jnr

As children grow and develop, a key feature of a children's wheelchair is that it can be adjusted to accommodate the child's individual growth and needs. With the Invacare Action<sup>3</sup> Junior, the width of the seat and backrest can be adjusted, as can the height of the arm and footrests.

A wide range of options are available on the Action<sup>3</sup> Junior to meet every child's needs in terms of comfort and positioning.

## Features and Options



### Seat-width adjustment

The seat width of the chassis can be adjusted with an innovative telescopic crossbar over 4 seat positions from 205 to 380 mm.




### Upholstery adjustment

As the chassis is made wider, so the upholstery on the back and on the seat of the chair can be made wider with Velcro flaps.



### Front stabilizer kit

A front stabilizer kit is available to increase forward stability for active children who are likely to climb all over the wheelchair and be especially energetic when propelling the wheelchair.



### Adjustable seat height and balance

For improved driving performance and different height options, the rear wheel can be mounted in 20 different positions.

## Tech Specs



**Seat width**  
20 - 38cm



**Seat depth**  
25 - 37.5cm



**Seat height**  
43.5cm (self-propel)  
46cm (transit)



**Backrest height**  
37, 40cm



**Total width**  
seat width + 17cm



**Total length**  
96cm



**Max. user weight**  
80kg



**Transport weight**  
8.5kg (self-propel)  
12kg (transit)



**Total weight unloaded**  
15kg



**Turning radius**  
75cm



**Max. safe slope**  
12°



FOR FURTHER INFORMATION:

**New Zealand:** Freephone 0800 468 222, email [sales@invacare.co.nz](mailto:sales@invacare.co.nz) or visit [www.invacare.co.nz](http://www.invacare.co.nz)

**Australia:** Freephone 1800 460 460, email [sales@invacare.com.au](mailto:sales@invacare.com.au) or visit [www.invacare.com.au](http://www.invacare.com.au)

**Asia:** Phone +61 2 8839 5330, email [exportsales@invacare.com.au](mailto:exportsales@invacare.com.au) or visit [www.invacare.com](http://www.invacare.com)

**Thailand:** Phone +66 0 2108 8555 or visit [www.invacare.com](http://www.invacare.com)



Yes, you can.®